



**ANALYSETTE 22**

**Prot.-No. M190301**

pre disp.in 1% Na4P2O7; 1min L-17; measured in 01% Na4P2O7

**Material CaCO3**




Q3(x) in %	x in µm	CV in %	M138679	M138699	M138700	M138723	M138724	M138744	M138745
10	0,98	127,6	4,01	0,69	0,69	0,39	0,39	0,34	0,34
50	6,11	176,1	32,3	3,59	3,54	0,9	0,91	0,76	0,76
90	29,91	186,8	165,65	18,37	17,98	1,99	1,99	1,7	1,7

x in µm	Q3(x) in %	CV in %	M138679	M138699	M138700	M138723	M138724	M138744	M138745
0,1	0,0	0,0	0	0	0	0,1	0,1	0	0
1,0	40,3	60,7	3,1	17,8	17,7	56,1	55,9	65,8	65,8
5,0	74,8	42,6	11,4	56,8	57,1	99	99	100	100
10,0	80,1	35,0	18,2	71,8	72,3	99,3	99,4	100	100
15,0	84,9	29,4	26	83,9	84,6	99,9	99,9	100	100
20,0	88,5	25,3	34,2	92,3	92,8	100	100	100	100
30,0	92,3	19,8	47,6	99,1	99,3	100	100	100	100
50,0	94,7	13,8	62,7	100	100	100	100	100	100
80,0	96,4	9,2	74,6	100	100	100	100	100	100
100,0	97,1	7,2	79,9	100	100	100	100	100	100
150,0	98,3	4,2	88,3	100	100	100	100	100	100
200,0	99,0	2,5	92,9	100	100	100	100	100	100
300,0	99,7	0,7	98	100	100	100	100	100	100
400,0	100,0	0,1	99,8	100	100	100	100	100	100
500,0	100,0	0,0	100	100	100	100	100	100	100

# ANALYSETTE 22

**Meas-No.** 138679 **Name** (3) **Date** 01-Okt-19 9:32:44  
**Material** CaCO3 **Lot** original sample  
**Information** Malaysia **Comment** M190301

**Calculation** Fraunhofer **TradeOff** broad (1000,0)  
**Refractive index** --- **Absorptions coefficient** ---  
**Scans Fine** 100 **Scans Coarse** 100 **Channels** 102 **Beam Obscuration** 13,0 %  
**Meas. Range** 0.08 - 2100 µm **Pump** 70 % **Ultrasonic** 100 %



138679 dQ3(x)
138679 Q3(x)


x in µm	Q3(x) in %
0,50	1,1
1,00	3,1
2,00	6,0
5,00	11,4
10,00	18,2
20,00	34,2
30,00	47,6
40,00	56,6
50,00	62,7
60,00	67,4
70,00	71,3
80,00	74,6
90,00	77,5
100,00	79,9

Q3(x) in %	x in µm
5,0	1,61
<b>10,0</b>	<b>4,01</b>
15,0	7,70
20,0	11,23
30,0	17,41
40,0	23,92
<b>50,0</b>	<b>32,30</b>
60,0	45,23
70,0	66,41
80,0	100,29
<b>90,0</b>	<b>165,65</b>
95,0	231,19
97,0	272,46
99,0	338,81

# ANALYSETTE 22

**Meas-No.** 138699 **Name** (2) **Date** 02-Okt-19 10:50:07  
**Material** CaCO<sub>3</sub> **Lot** 10min P5pl 20mm ZrO<sub>2</sub> in IPA  
**Information** Malaysia **Comment** M190301

**Calculation** Fraunhofer **TradeOff** broad (1000,0)  
**Refractive index** --- **Absorptions coefficient** ---  
**Scans Fine** 100 **Scans Coarse** 100 **Channels** 102 **Beam Obscuration** 12,0 %  
**Meas. Range** 0.08 - 2100 µm **Pump** 70 % **Ultrasonic** 100 %



138699 dQ3(x)
 — 138699 Q3(x)


x in µm	Q3(x) in %
0,10	0,0
0,50	5,2
1,00	17,8
2,00	35,5
3,00	45,8
4,00	52,4
6,00	60,4
8,00	66,3
10,00	71,8
12,00	77,0
15,00	83,9
20,00	92,3
30,00	99,1
40,00	100,0

Q3(x) in %	x in µm
5,0	0,49
<b>10,0</b>	<b>0,69</b>
15,0	0,88
20,0	1,10
30,0	1,63
40,0	2,38
<b>50,0</b>	<b>3,59</b>
60,0	5,89
70,0	9,34
80,0	13,25
<b>90,0</b>	<b>18,37</b>
95,0	22,45
97,0	25,07
99,0	29,65

# ANALYSETTE 22

**Meas-No.** 138700 **Name** (3) **Date** 02-Okt-19 10:50:40  
**Material** CaCO3 **Lot** 10min P5pl 20mm ZrO2 in IPA  
**Information** Malaysia **Comment** M190301

**Calculation** Fraunhofer **TradeOff** broad (1000,0)  
**Refractive index** --- **Absorptions coefficient** ---  
**Scans Fine** 100 **Scans Coarse** 100 **Channels** 102 **Beam Obscuration** 12,0 %  
**Meas. Range** 0.08 - 2100 µm **Pump** 70 % **Ultrasonic** 100 %



138700 dQ3(x)
  138700 Q3(x)


x in µm	Q3(x) in %
0,10	0,0
0,50	5,2
1,00	17,7
2,00	35,6
3,00	46,1
4,00	52,6
6,00	60,7
8,00	66,7
10,00	72,3
12,00	77,6
15,00	84,6
20,00	92,8
30,00	99,3
40,00	100,0

Q3(x) in %	x in µm
5,0	0,49
<b>10,0</b>	<b>0,69</b>
15,0	0,89
20,0	1,10
30,0	1,62
40,0	2,36
<b>50,0</b>	<b>3,54</b>
60,0	5,80
70,0	9,17
80,0	12,97
<b>90,0</b>	<b>17,98</b>
95,0	22,01
97,0	24,53
99,0	29,14

# ANALYSETTE 22

**Meas-No.** 138723 **Name** (3) **Date** 07-Okt-19 11:13:03  
**Material** CaCO<sub>3</sub> **Lot** 30min P5pl 0,5mm ZrO<sub>2</sub> in IPA  
**Information** Malaysia **Comment** M190301

**Calculation** Fraunhofer **TradeOff** broad (1000,0)  
**Refractive index** --- **Absorptions coefficient** ---  
**Scans Fine** 100 **Scans Coarse** 100 **Channels** 102 **Beam Obscuration** 17,0 %  
**Meas. Range** 0.08 - 2100 µm **Pump** 70 % **Ultrasonic** 100 %



138723 dQ3(x)
138723 Q3(x)


x in µm	Q3(x) in %
0,10	0,1
0,30	4,2
0,50	18,7
1,00	56,1
1,50	78,6
2,00	90,2
3,00	97,8
5,00	99,0
8,00	99,1
10,00	99,3
15,00	99,9
20,00	100,0
25,00	100,0
30,00	100,0

Q3(x) in %	x in µm
5,0	0,31
<b>10,0</b>	<b>0,39</b>
15,0	0,45
20,0	0,52
30,0	0,64
40,0	0,76
<b>50,0</b>	<b>0,90</b>
60,0	1,07
70,0	1,27
80,0	1,55
<b>90,0</b>	<b>1,99</b>
95,0	2,42
97,0	2,77
99,0	4,26

# ANALYSETTE 22

**Meas-No.** 138724 **Name** (4) **Date** 07-Okt-19 11:13:36  
**Material** CaCO<sub>3</sub> **Lot** 30min P5pl 0,5mm ZrO<sub>2</sub> in IPA  
**Information** Malaysia **Comment** M190301

**Calculation** Fraunhofer **TradeOff** broad (1000,0)  
**Refractive index** --- **Absorptions coefficient** ---  
**Scans Fine** 100 **Scans Coarse** 100 **Channels** 102 **Beam Obscuration** 16,0 %  
**Meas. Range** 0.08 - 2100 µm **Pump** 70 % **Ultrasonic** 100 %



x in µm	Q3(x) in %
0,10	0,1
0,30	4,1
0,50	18,6
1,00	55,9
1,50	78,5
2,00	90,1
3,00	97,8
5,00	99,0
8,00	99,1
10,00	99,4
15,00	99,9
20,00	100,0
25,00	100,0
30,00	100,0


Q3(x) in %	x in µm
5,0	0,32
<b>10,0</b>	<b>0,39</b>
15,0	0,46
20,0	0,52
30,0	0,64
40,0	0,77
<b>50,0</b>	<b>0,91</b>
60,0	1,07
70,0	1,28
80,0	1,55
<b>90,0</b>	<b>1,99</b>
95,0	2,43
97,0	2,77
99,0	6,00



# ANALYSETTE 22

**Meas-No.** 138744 **Name** (1) **Date** 08-Okt-19 13:31:48  
**Material** CaCO3 **Lot** 45min 800rpm P5pl 1mm ZrO2 in IPA  
**Information** Malaysia **Comment** M190301

**Calculation** Fraunhofer **TradeOff** broad (100,0)  
**Refractive index** --- **Absorptions coefficient** ---  
**Scans Fine** 100 **Scans Coarse** 0 **Channels** 51 **Beam Obscuration** 13,0 %  
**Meas. Range** 0.08 - 42 µm **Pump** 40 % **Ultrasonic** 100 %



138744 dQ3(x)
  138744 Q3(x)


x in µm	Q3(x) in %
0,10	0,0
0,20	0,8
0,30	6,7
0,40	16,0
0,50	26,2
0,60	36,0
0,80	52,9
1,00	65,8
1,20	75,6
1,50	85,5
1,80	91,6
2,00	94,3
3,00	99,5
4,00	100,0
5,00	100,0

Q3(x) in %	x in µm
5,0	0,28
<b>10,0</b>	<b>0,34</b>
15,0	0,39
20,0	0,44
30,0	0,54
40,0	0,64
<b>50,0</b>	<b>0,76</b>
60,0	0,90
70,0	1,08
80,0	1,32
<b>90,0</b>	<b>1,70</b>
95,0	2,07
97,0	2,31
99,0	2,76

# ANALYSETTE 22

**Meas-No.** 138745 **Name** (2) **Date** 08-Okt-19 13:32:03  
**Material** CaCO3 **Lot** 45min 800rpm P5pl 1mm ZrO2 in IPA  
**Information** Malaysia **Comment** M190301

**Calculation** Fraunhofer **TradeOff** broad (100,0)  
**Refractive index** --- **Absorptions coefficient** ---  
**Scans Fine** 100 **Scans Coarse** 0 **Channels** 51 **Beam Obscuration** 13,0 %  
**Meas. Range** 0.08 - 42 µm **Pump** 40 % **Ultrasonic** 100 %



138745 dQ3(x)
— 138745 Q3(x)

x in µm	Q3(x) in %
0,20	0,8
0,30	6,7
0,40	16,0
0,50	26,2
0,60	36,0
0,80	52,8
1,00	65,8
1,20	75,6
1,50	85,5
1,80	91,6
2,00	94,4
3,00	99,5
4,00	100,0
5,00	100,0

Q3(x) in %	x in µm
5,0	0,28
<b>10,0</b>	<b>0,34</b>
15,0	0,39
20,0	0,44
30,0	0,54
40,0	0,64
<b>50,0</b>	<b>0,76</b>
60,0	0,90
70,0	1,08
80,0	1,32
<b>90,0</b>	<b>1,70</b>
95,0	2,06
97,0	2,31
99,0	2,76