

Meas. No. 974 SOP 33 Date 04.08.2008 13:43:40 Operator FRITSCHLAN\Gerber

Material: **CAL Ca(OH)₂** in water

Description: in water + Na₄P₂O₇, 30s ultrasonic

Instrument parameter

Mode	Wet	Serial No.	22.3000/00835
Beam absorption	7.4 %	Stirrer	1.00 %
Meas. range	0.30 µm - 300.00 µm	Ultrasonic	On
Channels	62		
Scans	4	Prot.-No.:	080247

5-99%

Obere Kornklasse	Q3(x) [%]
1.033	5.0
1.558	10.0
2.061	15.0
2.572	20.0
3.103	25.0
3.657	30.0
4.240	35.0
4.856	40.0
5.514	45.0
6.228	50.0
7.013	55.0
7.886	60.0
8.875	65.0
10.017	70.0
11.371	75.0
13.026	80.0
15.137	85.0
18.189	90.0
23.223	95.0
34.131	99.0

0.3-300µm

Obere Kornklasse	Q3(x) [%]
0.330	0.3
0.500	1.1
1.000	4.7
2.000	14.4
3.000	24.0
5.000	41.1
8.000	60.6
10.000	69.9
15.000	84.7
20.000	92.2
30.000	98.1
40.000	99.7
50.000	100.0
60.000	100.0
80.000	100.0
100.000	100.0
150.000	100.0
200.000	100.0
250.000	100.0
300.000	100.0

