

Mess-Nr. 1037 SOP 33 Datum 17.09.2008 10:01:39 Benutzer FRITSCHLAN\Gerber

Material: **HP1 90µm**

Nassmessung

Dispergierung: in water + 0,1% Na4P2O7, 30s ultrasonic

Geräteparameter

Methode Wet

Serien Nr. 22.3000/00835

Strahlabsorption 12.6 %

Rührer 1

Meßbereich 0.30 µm - 300.00 µm

Ultraschall An

Kanäle 62

Scans 4

Prot.-Nr.: **080276**

5-99%

Obere Kornklasse [µ]	Q3(x) [%]
0.936	5.0
1.384	10.0
1.899	15.0
2.545	20.0
3.376	25.0
4.438	30.0
5.912	35.0
8.082	40.0
11.108	45.0
14.540	50.0
17.961	55.0
21.329	60.0
24.733	65.0
28.295	70.0
32.167	75.0
36.567	80.0
42.063	85.0
49.405	90.0
61.525	95.0
87.910	99.0

Kundenspezifische µm-Wert

Obere Kornklasse [µ]	Q3(x) [%]
30.000	72.3
63.000	95.5
71.000	97.2
90.000	99.1
100.000	99.6
150.000	100.0

