

Mess-Nr. 1075 SOP 34 Datum 14.10.2008 13:37:48 Benutzer FRITSCHLAN\Gerber

Material:

Blei/Zink-Pulver

Trockenmessung, 4bar

Geräteparameter

Methode	Dry	Serien Nr.	22.3000/00835
Strahlabsorption	2.7 %	Amplitude	0.00 %
Meßbereich	2.00 µm - 300.00 µm	Menge	0.00 %
Kanäle	31	Prot.-Nr.:	080301
Scans	3		

5-99%

Size Class	Q3(x) [%]
2.443	5.000
2.905	10.000
4.202	15.000
7.736	20.000
10.899	25.000
14.386	30.000
18.470	35.000
23.286	40.000
29.230	45.000
35.847	50.000
43.188	55.000
52.093	60.000
62.481	65.000
74.249	70.000
87.326	75.000
101.854	80.000
118.452	85.000
139.281	90.000
179.088	95.000
274.453	99.000

2-300µm

Size Class	Q3(x) [%]
2.000	5.151
5.000	16.163
10.000	23.621
15.000	30.784
20.000	36.716
30.000	45.594
40.000	52.882
50.000	58.871
60.000	63.891
80.000	72.264
100.000	79.421
125.000	86.712
150.000	91.727
200.000	96.168
250.000	98.070
300.000	99.972

