

Mess Nr. 5470 SOP 2 Datum 16.12.2008 08:04:15 Benutzer FRITSCHLAN\benes

Material: **Milchpulver**

dry measurement, 4bar

Dispergierung: pressure: 4 bar, concentration 1-4%

Methode Dry

Serien Nr. 22.2000.00/90771

Strahlabsorption 3,1 %

Zellpositionen 1

Meßbereich 4,45 µm - 610,69 µm

Berechnung Automatische Modellerkennung

Kanäle 51

Scans 100

Prot.-Nr.: **080364**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
31,078	5,0
46,534	10,0
56,520	15,0
64,980	20,0
72,457	25,0
79,601	30,0
86,612	35,0
93,586	40,0
100,596	45,0
107,687	50,0
115,315	55,0
123,308	60,0
131,697	65,0
140,992	70,0
151,497	75,0
163,426	80,0
177,630	85,0
196,441	90,0
224,809	95,0
278,020	99,0

5-600µm

Obere Kornklasse [µm]	Q3(x) [%]
5,000	0,4
10,000	2,4
15,000	3,1
20,000	3,5
30,000	4,8
50,000	11,6
60,000	17,0
80,000	30,3
100,000	44,6
120,000	58,0
150,000	74,3
200,000	90,7
250,000	97,5
300,000	99,5
350,000	99,9
400,000	100,0
450,000	100,0
500,000	100,0
550,000	100,0
600,000	100,0

