

Mess Nr. 5471 SOP 2 Datum 16.12.2008 08:11:01 Benutzer FRITSCHLAN\benes

Material: **Milchpulver**

dry measurement, 4bar

Dispergierung: pressure: 4 bar, concentration 1-4%

Methode Dry

Serien Nr. 22.2000.00/90771

Strahlabsorption 3,6 %

Zellpositionen 1

Meßbereich 4,45 µm - 610,69 µm

Berechnung Automatische Modellerkennung

Kanäle 51

Scans 100

Prot.-Nr.: **080364**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
31,400	5,0
46,711	10,0
56,763	15,0
65,310	20,0
72,867	25,0
80,083	30,0
87,169	35,0
94,219	40,0
101,308	45,0
108,517	50,0
116,177	55,0
124,257	60,0
132,773	65,0
142,055	70,0
152,695	75,0
164,784	80,0
179,189	85,0
198,255	90,0
226,402	95,0
279,749	99,0

5-600µm

Obere Kornklasse [µm]	Q3(x) [%]
5,000	0,4
10,000	2,3
15,000	2,9
20,000	3,3
30,000	4,7
50,000	11,5
60,000	16,8
80,000	29,9
100,000	44,1
120,000	57,4
150,000	73,8
200,000	90,4
250,000	97,4
300,000	99,5
350,000	99,9
400,000	100,0
450,000	100,0
500,000	100,0
550,000	100,0
600,000	100,0

