

Meas. No. 7193 SOP 0 Date 29.06.2009 16:14:25 Operator: FRITSCHLAN\Benes

Material: **LiCoPO4 (4)**

P7pl ZrO2 60min (Mono89)

Description: 2min ultrasonic + water + 0,1%Na4P2O7

Calculation Monomodal

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 1,2 %

Pump 65,00 %

Meas. range 0,01 µm - 53,10 µm

Cellposition 1

Channels 108

Ultrasonic On

Scans 100

Prot.-No.: **090124**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,011	5,0
0,013	10,0
0,014	15,0
0,016	20,0
0,018	25,0
0,020	30,0
0,022	35,0
0,024	40,0
0,026	45,0
0,028	50,0
0,030	55,0
0,032	60,0
0,035	65,0
0,038	70,0
0,042	75,0
0,046	80,0
0,051	85,0
0,058	90,0
0,070	95,0
0,101	99,0

001-1µm

Obere Kornklasse [µm]	Q3(x) [%]
0,010	4,6
0,015	16,8
0,020	30,9
0,025	43,6
0,030	55,1
0,035	64,7
0,040	72,5
0,045	78,9
0,050	83,9
0,055	87,9
0,060	91,1
0,070	95,0
0,080	97,2
0,090	98,4
0,100	98,9
0,500	99,9
1,000	100,0
1,500	100,0
2,000	100,0

