

Meas. No. 7358    SOP 45    Date 11.08.2009 11:46:07    Operator: FRITSCHLAN\Benes

Material: **Leaded glass powder**

**P7pl ZrO2 1mm 30min**

Description: 2min ultrasonic + water + 0,1%Na4P2O7

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 13,7 %

Pump 75,00 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **090158**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,249	5,0
0,324	10,0
0,395	15,0
0,466	20,0
0,540	25,0
0,620	30,0
0,706	35,0
0,802	40,0
0,911	45,0
1,036	50,0
1,183	55,0
1,358	60,0
1,576	65,0
1,851	70,0
2,214	75,0
2,749	80,0
3,540	85,0
4,915	90,0
7,392	95,0
12,146	99,0

01-30µm

Obere Kornklasse [µm]	Q3(x) [%]
0,120	0,1
0,200	2,3
0,300	8,3
0,500	22,3
1,000	48,6
1,500	63,5
2,000	72,2
3,000	81,8
4,000	87,0
5,000	90,2
6,000	92,6
7,000	94,4
8,000	95,8
9,000	96,9
10,000	97,8
15,000	99,7
20,000	100,0
25,000	100,0
30,000	100,0

