

Meas. No. 7868 SOP 36 Date 16.10.2009 09:42:37 Operator FRITSCHLAN\Benes

Material: BT

P7pl ZrO2 2mm 60min

Description: dispersed in isopropyl alcohol (2min ultrasonic)

Calculation Automatische Modellerkennung

Mode SVDU

Serial No. 22.2000.00/90771

Beam absorption 14,2 %

Pump ca. 60 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic Off

Scans 100

Prot.-No.: 090243

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,120	5,0
0,142	10,0
0,168	15,0
0,198	20,0
0,231	25,0
0,272	30,0
0,317	35,0
0,371	40,0
0,433	45,0
0,509	50,0
0,598	55,0
0,705	60,0
0,841	65,0
1,013	70,0
1,242	75,0
1,556	80,0
2,024	85,0
2,806	90,0
4,387	95,0
8,248	99,0

01-20µm

Obere Kornklasse [µm]	Q3(x) [%]
0,100	3,4
0,150	11,6
0,200	20,3
0,300	33,3
0,400	42,4
0,500	49,5
0,600	55,1
0,700	59,8
0,800	63,6
1,000	69,7
1,200	74,2
1,500	79,3
2,000	84,8
2,500	88,3
3,000	90,9
4,000	94,2
5,000	96,1
6,000	97,4
8,000	98,9
10,000	99,6
15,000	100,0
20,000	100,0

