

Meas. No. 8078 SOP 0 Date 03.11.2009 09:20:41 Operator: FRITSCHLAN\Benes

Material: MnO₂ -MIE

P7pl 80ml ZrO₂ 0,5mm 60min

Description: water + 0,1% Na₄P₂O₇, 30s ultrasonic

Calculation Sehr Breit

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 12,6 %

Pump 60,00 %

Meas. range 0,01 µm - 53,10 µm

Cellposition 1

Channels 108

Ultrasonic On

Scans 100

Prot.-No.: **090252**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,089	5,0
0,104	10,0
0,116	15,0
0,127	20,0
0,138	25,0
0,148	30,0
0,158	35,0
0,169	40,0
0,180	45,0
0,191	50,0
0,203	55,0
0,217	60,0
0,232	65,0
0,248	70,0
0,268	75,0
0,292	80,0
0,322	85,0
0,364	90,0
0,431	95,0
0,579	99,0

001-2µm

Obere Kornklasse [µm]	Q3(x) [%]
0,010	0,0
0,015	0,0
0,020	0,0
0,025	0,1
0,030	0,3
0,035	0,5
0,040	0,6
0,045	0,7
0,050	0,8
0,055	0,8
0,060	0,9
0,070	1,5
0,080	2,9
0,090	5,3
0,100	8,4
0,500	97,6
1,000	100,0
1,500	100,0
2,000	100,0

