

Mess Nr. 7588 SOP 36 Datum 9/15/2009 2:02:01 PM Benutzer FRITSCHLAN\Benes

Material: **Lactose P14; 0,2mm**

mit SVDU in Benzin

Dispergierung:

Geräteparameter

Methode	SVDU	Serien Nr.	22.2000.00/90771
Strahlabsorption	17.1 %	Pumpe	ca. 60%
Meßbereich	0.77 µm - 504.48 µm	Zellpositionen	2
Kanäle	102	Ultraschall	Aus
Scans	100	Prot.-Nr.:	090254

5-99%

Size Class	Q3(x) [%]
1.843	5.000
3.772	10.000
5.980	15.000
7.971	20.000
9.783	25.000
11.507	30.000
13.212	35.000
14.907	40.000
16.637	45.000
18.445	50.000
20.340	55.000
22.382	60.000
24.614	65.000
27.329	70.000
30.340	75.000
34.155	80.000
39.478	85.000
48.313	90.000
71.818	95.000
117.419	99.000

01-250µm

Size Class	Q3(x) [%]
0.150	0.334
0.500	0.334
1.000	1.376
3.000	8.221
5.000	12.724
8.000	20.077
10.000	25.621
15.000	40.273
20.000	54.105
30.000	74.452
40.000	85.387
50.000	90.596
60.000	93.157
80.000	95.984
100.000	97.900
150.000	99.865
200.000	100.000
250.000	100.000

