

Meas. No. 8159 SOP 18 Date 16.11.2009 13:24:10 Operator: FRITSCHLAN\Benes

Material: **Microsilica**

P7pl WC 80ml 0,6mm 60min

Description: in H₂O + 0,1% Na₄P₂O₇ (30s Ultrasonic)

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 14,3 %

Pump 60,00 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **090271**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,182	5,0
0,243	10,0
0,305	15,0
0,375	20,0
0,450	25,0
0,536	30,0
0,632	35,0
0,744	40,0
0,876	45,0
1,029	50,0
1,215	55,0
1,439	60,0
1,714	65,0
2,062	70,0
2,522	75,0
3,108	80,0
3,904	85,0
5,072	90,0
7,054	95,0
10,978	99,0

01-30µm

Obere Kornklasse [µm]	Q3(x) [%]
0,120	0,7
0,200	6,4
0,300	14,6
0,500	28,0
1,000	49,1
1,500	61,3
2,000	69,2
3,000	79,1
4,000	85,5
5,000	89,7
6,000	92,7
7,000	94,9
8,000	96,5
9,000	97,6
10,000	98,4
15,000	99,9
20,000	100,0
25,000	100,0
30,000	100,0

