

**Mess.-Nr.** 60355 **Datum** Mittwoch, 30 Januar 2013 02:58 **Bearbeitungs- Nr.:**  
**Material** NaCl (250g - trocken) **P5cl 500ml 20mm ZrO2 20min** **M130007**  
**Bemerkung** predis. in IPA 1 min Ultrasonic

**Berechnung** Fraunhofer broad **TradeOff** 1.000000e+003  
**Scans** 100 **Scans** 100 **Messkanäle** 108 **Strahlabsorption** 17.00 %  
**Feinbereich** **Grobbereich**  
**Messbereich** 0,01 - 2000.00 [µm] **Pumpe** 5 **Ultraschall** 0

| %  | < µm  |
|----|-------|
| 5  | 0.86  |
| 10 | 1.20  |
| 15 | 1.55  |
| 20 | 1.94  |
| 25 | 2.41  |
| 30 | 3.01  |
| 40 | 4.88  |
| 50 | 8.15  |
| 60 | 12.50 |
| 70 | 17.71 |
| 80 | 24.23 |
| 90 | 33.69 |
| 95 | 41.49 |
| 98 | 50.29 |
| 99 | 55.88 |

| %     | < µm |
|-------|------|
| 0.1   | 0.1  |
| 1.0   | 0.5  |
| 7.0   | 1.0  |
| 40.5  | 5.0  |
| 54.5  | 10.0 |
| 65.0  | 15.0 |
| 73.8  | 20.0 |
| 81.0  | 25.0 |
| 86.6  | 30.0 |
| 91.0  | 35.0 |
| 94.2  | 40.0 |
| 97.9  | 50.0 |
| 99.4  | 60.0 |
| 99.9  | 70.0 |
| 100.0 | 80.0 |

