


Mess.-Nr. 60296 **Datum** Mittwoch, 30 Januar 2013 09:31 **Bearbeitungs- Nr.:**
Material NaCl (in 150ml Salatöl) **PI 5cl 400rpm 500ml ZrO2 30mm** **M130007**
Bemerkung predis. in IPA 1 min Ultrasonic

Berechnung Fraunhofer broad **TradeOff** 1.000000e+003
Scans 100 **Scans** 100 **Messkanäle** 108 **Strahlabsorption** 13.00 %
Feinbereich **Grobbereich**
Messbereich 0,01 - 2000.00 [µm] **Pumpe** 5 **Ultraschall** 0

%	< µm
5	0.69
10	0.93
15	1.14
20	1.37
25	1.61
30	1.88
40	2.56
50	3.64
60	5.72
70	9.36
80	14.09
90	20.82
95	26.57
98	33.00
99	37.07

%	< µm
0.3	0.1
2.2	0.5
11.7	1.0
32.1	2.0
57.3	5.0
71.5	10.0
81.6	15.0
89.0	20.0
93.9	25.0
96.9	30.0
98.6	35.0
99.4	40.0
99.8	45.0
100.0	50.0
100.0	60.0



 60296 dQ3(x) — 60296 Q3(x)