

Fritsch GmbH
Industriestrasse 8
55743 Idar-Oberstein
Germany



www.fritsch.de info@fritsch.de
Tel: 0049 6784 700 / Fax: 0049 6784 7011

**Fibratory Feeder LABORETTE 24 / Vibrations-Zuteilrinne LABORETTE 24 &
Rotary Cone Divider LABORETTE 27 / Rotations-Kegelprobenteiler LABORETTE 27**

Head / Kopf: 1:10
Sample / Probe: **Catalyst M130008** **Result 1-2**
Time (min): 15 fed centric to the cone

Bottle (No.)	Amount (g)
1,00	24,39
2,00	24,79
3,00	25,04
4,00	25,04
5,00	24,64
6,00	24,86
7,00	24,51
8,00	25,10
9,00	24,49
10,00	24,97

Average weight
Durchschnittsgewicht (g):
24,78

Standard deviation +/-
Standardabweichung +/-
0,247

total: **247,83**

Feed rate
Förderrate (g/min):
16,52