

Meas.Nr. 65413 **Date** Mittwoch, 26 Juni 2013 04:31 **Personal Report No.:**
Material Limestone 10min P6cl 30 mm ZrO2 **M130172**
Description predisp.. in 0,1% Na4P2O7 + 1 min LABORETTE 17 ultrasonic

Calculation Fraunhofer broad **TradeOff** 1.000000e+003
Scans Fine 100 **Scans Coarse** 100 **Channels** 102 **Beam Obscuration** 16.00 %
Meas. Range 0.08 [μm] - 2000.00 [μm] **Pump** 6 **Ultrasonics** 10

%	< μm
5	0.49
10	0.70
15	0.93
20	1.20
25	1.51
30	1.90
40	3.02
50	4.87
60	7.43
70	10.54
80	14.61
90	20.98
95	26.95
98	33.87
99	38.39

%	< μm
5.2	0.5
16.3	1.0
31.1	2.0
50.6	5.0
68.4	10.0
80.8	15.0
88.8	20.0
93.7	25.0
96.6	30.0
98.3	35.0
99.2	40.0
99.7	45.0
99.9	50.0
100.0	60.0

