

Meas.Nr.	65484	Date	06.28.2013 08:40			Personal Report No.:
Material	Limestone		1 h P6cl 5mm ZrO2			M130172
Description	predisp. in 0,1% Na4P2O7 with 1min ultrasound LABORETTE 17					
Calculation	Mie	n = 1.51000	-i 0	broad TradeOff	1.000000e+002	
Scans Fine	100	Scans Coarse	0	Channels	57	Beam Obscuration 13.00 %
Meas. Range	0.08 [µm]	- 42.30 [µm]	Pump	6	Ultrasonics	10

%	< µm
5	0.097
10	0.113
15	0.131
20	0.150
25	0.171
30	0.195
40	0.256
50	0.347
60	0.513
70	0.928
80	1.570
90	2.301
95	2.873
98	3.513
99	3.950

%	< µm
2.71	0.01
2.71	0.05
6.11	0.10
30.90	0.20
45.52	0.30
54.04	0.40
59.44	0.50
63.07	0.60
67.79	0.80
71.16	1.00
86.30	2.00
95.85	3.00
99.81	5.00
99.90	8.00
99.90	10.00

