


Mess.Nr. 69563 **Datum** 10/29/2013 15:20 **Bearbeitungs-Nr:**
Material BoSi-Glas (MTP-Messung) Mie Neuberechn. 90min P7pl 0,5 mm ZrO2 **M130270**
Bemerkung in Wasser + 0,1% NA4P2O7, ohne US

Berechnung Mie n = 1.50187 -i 0 **TradeOff** broad (100.0)
Scans 100 **Scans** **Messkanäle** 57 **Strahlabsorption** 11.00 %
Feinbereich **Grobbereich** 0
Messbereich 0.08 [µm] - 42.30 [µm] **Pumpe** 5 **Ultraschall** 10

%	<µm
5	0.089
10	0.090
15	0.096
20	0.101
25	0.108
30	0.115
40	0.130
50	0.149
60	0.173
70	0.205
80	0.252
90	0.334
95	0.429
98	0.588
99	0.810

%	<µm
9.5	0.08
18.9	0.10
50.3	0.15
68.6	0.20
86.8	0.30
93.8	0.40
96.8	0.50
98.1	0.60
98.7	0.70
99.0	0.80
99.2	0.90
99.3	1.00
99.6	1.20
99.8	1.50
100.0	2.00



▲ 69563 dQ3(x) — 69563 Q3(x)